



Handicare 2000

A Minimum swivel radius	690mm
– with arm spacers	713mm
– with extended arm spacers	730mm
B Minimum wall to seat back	95mm
C Seat depth	360mm
– with extended seat	425mm
D Minimum folded width to footplate	400mm
E Minimum open width to edge of footplate	689mm
F Armrest width – external	550mm
– with arm spacers	590mm
– with extended arm spacers	637mm
G Armrest width – internal	460mm
– with arm spacers	516mm
– with extended arm spacers	562mm
H Seat back height	453mm
I Footplate to seat height	
– standard seat option	485-630mm
– perch seat option	650-875mm
– fixed seat option	426-446mm
J Minimum footplate height	70mm
Bottom parking space required	700mm†
Minimum track intrusion into staircase	116mm
Minimum track intrusion at bottom of stairs	500mm
Minimum track intrusion at bottom of stairs – reduced lower intrusion	420mm
Maximum weight limit	115kg/18 stone
– heavy duty version	137kg/21.5 stone
Speed	0.12m/sec (approx)
Power supply	24v DC battery
Charger	100-240v (50-60Hz)
Maximum staircase angle	56°

† +20mm with arm spacers, +44mm with extended arm spacers.



User support: 08700 118282

Handicare, 82 First Avenue, Pensnett Estate, Kingswinford, West Midlands DY6 7FJ, United Kingdom
Textphone: (180015) 020 7490 8913 ext 05 Voicephone: (18002) 020 7490 8913 ext 05
Global Dealer Support: +44 (0) 844 225 3121 enquiries@handicare.co.uk www.handicare.co.uk

Please be aware that actual colours may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.



handicare

Handicare 2000

Curved stairlift



The benefits of a Handicare 2000 curved stairlift

If you are thinking about installing a stairlift, and your stairs turn corners or are spiral, the Handicare 2000 offers you a safe and reliable solution.

Handicare stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Handicare 2000 stairlift suits most staircases with bends or corners and offers you all the basic functions that you would expect plus the option to upgrade to a powered seat swivel and/or folding footplate.

With a choice of six upholstery colours and two track colours the Handicare 2000 compliments most colour schemes.

Fitted to the stairs, not the wall, it is unlikely you will need to redecorate following the installation of your Handicare 2000 stairlift.

Like all Handicare stairlifts the 2000 Series has a battery backup which means it will continue to work in the event of a power cut so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and the bespoke track is designed to fit your individual staircase and maximise clear space on the stairway. It can be supplied with a powered hinged track if required.

The Handicare 2000 also comes with two remote controls which means, once folded, you can easily send it to the upstairs charge point freeing up space in your hallway.



Safety as standard

Your Handicare 2000 stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

Your Handicare 2000 stairlift has an easy release retractable seat belt for added peace of mind. If you prefer you can have an easy fasten Velcro seat belt.

Easy to use

The operating toggle is easy to use, simply hold it towards the direction you wish to travel, up or down, and your Handicare 2000 will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.

